

MediCaption Readme

Github link: <https://github.com/aarthi03/MediCaption.git>

MediCaption Prototype Explanation:

We have created a sample UI to upload a Medicine Image for the medicine name detection.

Sample UI: For the sample UI we have developed a web page to upload images by user and storing it into a database.

We can access the UI from the link: <http://localhost:3000/>

We have uploaded an image:


Image Uploader



Choose File C-Retard 7.jpg

Upload Image

Image Uploader



Choose File C-Retard 7.jpg

Upload Image

Image uploaded successfully

Input image:



Initialization of server.js file which helps us connect html file with the database. For every file we upload we can see the DB updation from server.js

```
C:\Windows\System32\cmd.e  X  +  v
Microsoft Windows [Version 10.0.22621.3155]
(c) Microsoft Corporation. All rights reserved.

C:\Users\saira\OneDrive\Desktop\website\Test>node server.js
Server is running on port 3000
Connected to the database.
Table created successfully.
Image uploaded successfully with rowid 7
Image uploaded successfully with rowid 8
```

Initially we trained the model using: “YoloV8_TRAIN,VALID,TEST.ipynb”
From this file we got the best.pt with the model weights.

We have extracted the content of the database file by uploading the database into the git folder and we have also extracted the git folder files into the finalprototype.ipynb

```
Downloading the data from the git repository where the database is stored after updation

[ ] !git clone https://github.com/aarthio3/MediCaption.git

Cloning into 'MediCaption'...
remote: Enumerating objects: 3092, done.
remote: Counting objects: 100% (111/111), done.
remote: Compressing objects: 100% (101/101), done.
remote: Total 3092 (delta 23), reused 72 (delta 7), pack-reused 2981
Receiving objects: 100% (3092/3092), 1.61 GiB | 44.29 MiB/s, done.
Resolving deltas: 100% (69/69), done.
Updating files: 100% (2958/2958), done.

Double-click (or enter) to edit
```

The file best.pt is obtained while training the data using the code "YoloV8_TRAIN,VALID,TEST.ipynb" and then the file is stored in the GIT repository so that we can use the best.pt file which contains the weights of the model to predict the bounding boxes for the images that are being uploaded


```
Predicting the bounding boxes for the medicine image uploaded in the new images folder where the model predicting the labels probability is equal to or greater than 60%

!yolo predict model=/content/Test/best.pt source='/content/new_images' save_crop conf=0.6

Ultralytics YOLOv8.1.24 Python-3.10.12 torch-2.1.0+cu121 CPU (AMD EPYC 7B12)
Model summary (fused): 168 layers, 11125971 parameters, 0 gradients, 28.4 GFLOPs

image 1/1 /content/new_images/Ciprofar14.jpg: 640x480 1 Medicine Name, 913.7ms
Speed: 13.3ms preprocess, 913.7ms inference, 32.2ms postprocess per image at shape (1, 3, 640, 480)
Results saved to /content/ultralytics/runs/detect/predict2
Learn more at https://docs.ultralytics.com/modes/predict
```

From the image uploaded we have extracted the text of the medicine name by using OCR and Trained YoloV8 model :



Based on the Medicine Name which we have obtained from the above steps We have generated the details about the medicine and the corresponding audio file.

